

# Traveling with an IV Infusion Pump

## Planning & Packing Recommendations

For people with an IV infusion pump, travel can often cause extra stress. But with a little preparation, traveling while on IV therapy can be a safe and manageable experience. To help patients with easing their travel worries, here are a few things we suggest you prepare before traveling.

### Before booking travel (or at least 72 hours in advance):

- Confirm with your or your child's care team if travel is safe.
- Discuss with your care team any activities you plan to participate in to confirm they are safe.
- Locate a hospital or care provider at your destination. It may be helpful to ask your care team if they have recommendations.
- Discuss with your care team what supplies they recommend you bring and the quantity of each supply. We have included space below to take any needed notes.
- If transporting parenteral nutrition (PN), determine which type of cooler you will need. Your provider may have recommendations.
- Contact TSA through the TSA Cares program to make arrangements and ask any questions that you may have. Be sure to do this at least 72 hours before your travel. Call (855) 787-2227 or visit <https://www.tsa.gov/contact-center/form/cares>
- Contact your airline to notify them that you will be carrying medically necessary supplies on-board.
- Contact your hotel or vacation rental to see if shipping supplies to your destination is possible and to ensure a refrigerator is available at your destination.
- Determine if you will be carrying all your nutrition supplies with you or if they can be shipped directly to your destination.
- Download the CURLIN Travel Letter from Moog's website: <https://www.moogmedical.com/air-travel-curlin-infusion-pump/>

### The day prior to travel:

- Pack your PN or IV fluids into coolers with the appropriate cooling packs, as directed by your care team.
- Pack any emergency supplies recommended by your care team.

### CURLIN pump accessories to bring:

If you are traveling with an CURLIN Infusion Pump, we recommend bringing the following accessories with you in your carry-on bag or suitcase:

- Infusion pump(s)
- Batteries for infusion pump
- AC adapter for infusion pump
- IV pole clamp holster and IV pole (or an IV pole alternative, such as FreeArm®) or a pump carry pack
- PN or IV hydration tubing
- Insurance information
- Contact information for your home care provider and care team.
- Extra PN or IV fluids, and any additional required medication (check with your care team to determine how many extra days of medication to bring).

### Notes: