

7/16/2020

Dear Customer,

Below you will find the product description provided to Moog Medical from the manufacturer of the 1.2 and 0.22 micron filters on all CURLIN filtered administration sets. This pertains only to the filter properties, and not the administration set as a whole.

Filter Specifications:

- Fluid filter material: Hydrophilic membrane, Polyethersulfone
 - Filter membrane is low protein binding
 - Filter membrane is surfactant free
- Air filter material: Hydrophobic membrane, Polytetrafluoroethylene
 - Pore size 0.02 micron
- Non-DEHP
- BPA free
- Material formulation not to include natural rubber latex

Please feel free to contact the Clinical Services Department should you require further information.

Regards,
Clinical Services
Moog Medical