



Device Names:

- Infinity Enteral Feeding Pump
- Infinity Orange Enteral Feeding Pump

Manufacturer: Moog, Inc.

An enteral feeding pump is a medical device used to deliver liquid nutrition directly into a patient’s stomach or intestines through a feeding tube. It ensures that the nutrition is administered at a controlled rate to help maintain proper nutrition and hydration for individuals who cannot eat by mouth due to various medical conditions.

- This device is medically necessary.
- This device complies with TSA security screening requirements for external medical devices.
- This device is non-transmitting and does not contain Bluetooth or WiFi technology.
- This device has been tested to the aviation standards RTCA DO/160E parts 20 and 21 for radiated RF energy and susceptibility and IEC 60601-1 for safety.

TSA Security Screening for External Medical Devices

It is recommended that the Moog Infinity and Moog Infinity Orange Feeding pump be screened using:

- Metal Detector
- Personal Physical Screening

*** Power off pump prior to passing it through the metal detector.**

The Infinity and Infinity Orange feeding pumps use a rechargeable lithium-ion battery, which complies with the FAA’s policy on lithium batteries.

Does device contain a lithium-ion battery?	Yes
Is lithium-ion battery spare or installed?	Installed
Battery complies with FAA watt hour (Wh) rating policy?	Yes
Moog Infinity battery watt hour (Wh) rating	9.62 watt hours (Wh)



bit.ly/3Yqy63q

Scan the code above to see the FAA’s policy on external medical devices.



bit.ly/4hCEEoy

Scan the code above to see the FAA’s policy on lithium batteries.