

NOTICE: New CURLIN® IV Administration Sets Arriving Soon

Dear Valued Customer,

This notice is to inform you of an update for CURLIN Administration Sets. Moog is discontinuing the 340-XXXX CURLIN Administration Sets (340 Series) * and launching a new family of CURLIN administration sets: the 380-XXXX CURLIN Administration Sets (380 Series). These new 380 administration sets will have all the same features and configurations as the current 340 administration sets, but with some enhancements (see "Enhancements" below for more detail).

**exception for 340-4115 and 340-4133, which will continue to be orderable*

Timing:

We will no longer accept orders for the 340 administration sets beginning on November 3rd, 2025. For planning purposes, please consider that while we will no longer accept new orders for the 340 administration sets, we will continue fulfilling orders for the 340 administration sets that were placed before November 3rd, 2025.

Enhancements:

- Backwards compatibility: The 380 administration sets will be available for use in both the CURLIN 8000 and CURLIN 6000.
- Easier set loading: the yellow Free-Flow Protection housing will have a longer arm to assist in loading the set into the pump.
- 96-hour set life: the 380 administration sets have been validated for up to 96 hours of use with the CURLIN 8000. *Note: when using the 380 administration sets with the CURLIN 6000, they are validated for 72 hours*

New Filters:

- The 380 filtered administration sets (both 0.2 and 1.2 micron) will have a different filter than the filters seen on the 340 administration sets. The function is the same, however, priming instructions are different. Please be sure to read the [priming instructions included in the packaging for each set](#) and educate your end users accordingly.

Ordering Process:

When the 380 administration sets are released into the market, the ordering process will remain the same as the current ordering process. Please continue to order your sets through your distributor or directly through Moog Medical Customer Service.

When placing an order, the prefix for each new administration set will be '380' (the prefix for the current administration sets is '340'). The suffix for each updated product code will remain the same as the suffix in the current product line. For example, the current product code 340-4114 will become 380-4114. Refer to the table below for a list of the new 380 product codes that are equivalent to 340 product codes or new codes altogether. Equivalent product codes are in **green** text, new codes are in **blue** text.

"340" Product Codes being discontinued		"380" Product Codes direct replacements
340-4114	▶	380-4114
340-4126	▶	380-4126
340-4128	▶	380-4128
340-4128-V	▶	380-4128-V
340-4130	▶	380-4130
340-4130-V	▶	380-4130-V
340-4134	▶	380-4134
340-4144	▶	380-4144
340-4165	▶	380-4165
340-4166	▶	380-4166
340-4168	▶	380-4168
340-4173	▶	380-4173
340-4174	▶	380-4174
New 380 Administration Set configurations		
		380-4178
		380-4192
		380-4194
		380-4195

Some of the 380 administration sets are completely new sets. Refer to the [380 set brochure](#) or our [website](#) for more details and a quick reference guide for the new product codes. Lastly, **please consider the following:**

- **Pricing for the new 380 administration sets will be consistent with your current pricing for equivalent 340 administration sets.**

- **Any orders placed for 340 administration sets after November 3rd, 2025, will automatically be fulfilled with the 380 administration set equivalent.**
- **To ensure sufficient supply to all customers during this transitional period, we will not honor any large stocking orders outside of typical order volumes for either 340 or 380 administration sets.**

The 380 administration sets will be backed by the same industry-leading clinical and customer support teams, which are available to assist customers and patients with any question or concern 24 hours a day, 7 days per week, any day of the year.

We truly appreciate your business. If you have any further questions related to the new CURLIN 380 Series Administration Sets, please contact your Moog Sales Representative.

Sincerely,

Alex Merrill

Sr. Product Manager, Infusion

A decorative graphic in the bottom left corner of the page, consisting of several overlapping, flowing, light blue and white shapes that create a sense of movement and depth.